




RaspberryMed

Hybrid Dilator Storage Cart



- 24"W x 40"D X 45"H
- Clear acrylic doors w/ lock
- 120 dilator capacity
- 6 trays hold up to 12 dilators each, for a total of 72 dilators stored **horizontally**
- 4 cassettes hold up to 12 dilators each, for a total of 48 dilators stored **vertically**
- Drawer
- 4 locking casters
- Available 402NM LED lighting
- Available HEPA filtering system
- Available keyless lock upgrade



RaspberryMed

For more information email us at sales@raspberrymed.com
Toll-free 800-213-9704 | Telephone (616) 949-0080 | Fax (616) 949-8333