




Endoscope Storage Cart



- 38.75"W x 77.75"T x 28"D
- Holds 8 scopes
- Hanging up to 63" standard
- Up to 77" hanging length with optional roller brackets
- FDA approved silicone padded back wall
- FDA approved nitrile soft grip scope holders
- No maintenance hinged acrylic doors
- Cam lock (keyed), Standard
- Electronic keyless lock (optional)
- 5" locking casters
- Lightweight aluminum frame
- Superior stability and ergonomic scalloped base
- HEPA filtering system, standard
- 405nm disinfecting light, safe use for endoscopes



RaspberryMed

For more information email us at sales@raspberrymed.com

Toll-free 800-213-9704 | Telephone (616) 949-0080 | Fax (616) 949-8333