

## WALL MOUNT VERTICAL DILATOR STORAGE:

19" W X 4" D X 36" T (24 Dilator Storage)




### Features


- Slim design allows for 24 or 48 vertical hung dilators depending on model
- Patent Pending, FDA compliant soft grip system for 5mm to 20mm dilators
- Key locking system with available keyless upgrade
- Available 405nm disinfecting light and/or HEPA filtering system



38" W X 4" D X 38" T (48 Dilator Storage)



FOR PRICING AND DETAILED INFORMATION CONTACT  
RASPBERRY MED / 616.949.0080  
[jeff@raspberrymed.com](mailto:jeff@raspberrymed.com)



RaspberryMed

For more information email us at [sales@raspberrymed.com](mailto:sales@raspberrymed.com)  
Toll-free 800-213-9704 | Telephone (616) 949-0080 | Fax (616) 949-8333


# Medley Dilator Storage

## **MOBILE VERTICAL DILATOR STORAGE:** 24" W X 24" D X 44" T (72 Dilator Storage)



### **Features**

- Patent pending, FDA compliant soft grip system for 5mm to 20mm dilators
- Key locking system with available keyless upgrade
- Available 405nm disinfecting light
- Available HEPA Filtering System
- Slide out mounting storage trays



RaspberryMed

For more information email us at [sales@raspberrymed.com](mailto:sales@raspberrymed.com)  
Toll-free 800-213-9704 | Telephone (616) 949-0080 | Fax (616) 949-8333