



RaspberryMed

KEY FEATURES

- Dimensions: 38.75" W x 28" D x 77.75" T
- Holds 8 Scopes
- Up to 77" hanging length with optional roller brackets
- Hanging up to 63" as standard.
- FDA approved silicone padded back wall
- FDA approved nitrile soft grip scope holders
- No maintenance hinged acrylic doors
- Key Lock or optional Electronic Keypad Lock
- 5" Locking Casters
- Lightweight all aluminum frame
- Superior stability and ergonomic scalloped base
- HEPA filter standard
- 405nm disinfecting lighting (Optional) Safe for use on Endoscopes



Optional Rollers Shown

