





RaspberryMed

Esophageal Dilator Storage Wall Mount



- Single (holds 24 dilators) - 19"W x 36"H x 4"D
- Double (holds 36 dilators) – 38"W x 38"H x 4"D
- Patent pending, FDA compliant soft grip attachment system for 5mm to 20mm dilators
- Cam lock, standard
- HEPA filtering system (optional)
- 405nm disinfecting light (optional)
- Electronic keyless lock (optional)



RaspberryMed

For more information email us at sales@raspberrymed.com

Toll-free 800-213-9704 | Telephone (616) 949-0080 | Fax (616) 949-8333