



RaspberryMed

## Respiratory Cart

26.75"W x 32.8"D x 55"H (plus IV pole height)



- Respiratory cart (ECMO / ECLS / ELSO) ideal for ventilators, oxygen tanks
- Stainless 304 construction
- 1.25"-diameter post for standard medical equipment attachments
- 5"-diameter locking casters with 2" wide wheels that can't get stuck in elevator threshold gaps
- Height-adjustable IV pole that can be positioned on any three corners
- Removable post for easy installation and removal of equipment oxygen tank holder, movable to any cart corner
- Six-plug patient vicinity grade outlet on shelf for easy cable management
- Upper and lower shelves can hold 90 lbs. each
- Recessed frame and outrigged casters provide low center of gravity for maximum stability
- Removable side handles



RaspberryMed

For more information email us at [sales@raspberrymed.com](mailto:sales@raspberrymed.com)

Toll-free 800-213-9704 | Telephone (616) 949-0080 | Fax (616) 949-8333