




RaspberryMed

## Enteroscope Storage Cart



- 24" W X 12"D x 102" H
- Fits All Enteroscopes: Hanging length up to 95"
- Stores 5 scopes per cabinet
- Pull down scope manager allows for easy and ergonomic access to scope by lowering the scopes 18 inches to a reachable height.
- HEPA filter and 405NM clean wave lighting keep in compliance with SGNA and other GI industry standards
- FDA Approved soft grip material and backing ensures the endoscope delicate ends and controls do not get damaged
- Easy two piece construction make the cabinet easy to transport and maneuver into place and makes installation quick and intuitive
- Full length acrylic door is maintenance free and allows quick visual inspection into the cabinet for inventory compliance
- Removable bottom tray for maximum cleanability
- Standard key lock with optional push button lock



RaspberryMed

For more information email us at [sales@raspberrymed.com](mailto:sales@raspberrymed.com)

Toll-free 800-213-9704 | Telephone (616) 949-0080 | Fax (616) 949-8333