



RaspberryMed

Point of Care (POC) Cart

22" W x 27.5" D X 38.5 "H



Features:

- Height adjustable work surface. Adjusts up/down 12" from 34" to 46" Tall
- Solid surface tops and worksurface areas for maximum cleanability
- Left and right slide out worksurface
- slide out keyboard tray with dual mouse pads
- Integrated drip edge on all top and slide out work surfaces
- Drawers: 12.75" W x 11.5" Deep in 3" and 6" tall versions
- Lockable laptop storage drawer
- Soft Close drawer slides
- Integrated push handle
- Three locking a one tracking caster. All 5" diameter
- Two lock zones available in key lock or electronic pushbutton lock
- Side rails for mounting accessories



RaspberryMed

For more information email us at sales@raspberrymed.com

Toll-free 800-213-9704 | Telephone (616) 949-0080 | Fax (616) 949-8333