



RASPBERRY MED



Point of Care (POC) Cart

PRODUCT OVERVIEW



Point of Care (POC) Cart



The Point of Care Cart brings charting, medications, and supplies to the bedside in one compact unit—height adjustable, lockable, and built with top tier quality and durability.

Features

Height-Adjustable Computer Station

The computer station adjusts a full 12 inches, from 34" to 46" in height, providing flexibility to accommodate staff of varying heights.

Power Options

Runs on standard AC power, with a rechargeable battery available as an option for units that stay on the move through a full shift.

360° Rotation

The computer station rotates a full 360°, locking every 90° so the screen can face the clinician.

Keyboard Tray

The tray slides out from under the worksurface and carries a mouse pad on each side, so the cart works for left and right handed staff.

Two Lock Zones

Drawer storage and the computer station secure independently, with a key lock or an electronic push-button.

10-Year Warranty

Covers materials and workmanship for a period of 10 years.



POINT OF CARE CART



Work Surface

The work surface slides left to right across the workstation to support both left and right-handed clinicians.



Mixed 3" and 6" Drawers

Specify any combination of 3" and 6" drawers, adding up to a total of 21" in height.



360° Rotation

The computer station rotates a full 360° and locks every 90° so the screen can face the clinician.



BUILT TO ORDER

Configure Your Cart

Most hospitals utilize two carts per room: one for charting, and one for supplies. The Raspberry Med Point of Care Cart covers both needs in one. With only a 20" x 25" footprint, it offers everything your care teams need in a compact, efficient design.

Drawers

The cart's smooth gliding drawers measure 16.5" by 11.5" and come in 3" and 6" heights. Customize a combination of drawer heights to your needs, adding up to a total of 21".

Locks

Two lock zones are built-in: one locks all of the drawers, the other locks the computer station. Lock options include a key or electronic push button.

Your Logo & Finishes

Select the right finish to match your facility's aesthetic, from our variety of paint colors and wood grain options. We can also print your logo on the side panels of the cart, bringing added brand recognition to every unit.


Side Rails

Optional accessory rails allow you to easily mount attachments to the cart, giving your care team portable, easy access to additional accessories and equipment.


Finishes

The Point of Care Cart comes in a variety of finishes and wood grain laminate options to match your hospital's aesthetic. Each finish is a medical-grade material, designed to be durable and easy to clean.

COLOR OPTIONS



White



Grey



Beige



Taupe

WOOD-GRAIN LAMINATE OPTIONS



Light Oak



Natural Walnut






Specs

Get all the functionality your team needs to care for patients in one all-in-one solution. The Raspberry Med Point of Care Cart is compact, adjustable, and easy to clean. Featuring gliding drawers with plenty of space, the POC is lockable and durable. Easily adjust the height, rotate the computer station a full 360 degrees, or slide the work surface left and right. The POC Cart is designed to adapt quickly to your needs, removing distractions so you can focus on patient care.

Overall dimensions	20.25"W x 25.07"D x 39.28"H
Worksurface height	34" to 46", 12" of travel
Worksurfaces	Left and right slide-out
Rotation	360°, locks every 90°
Keyboard tray	Slide-out with dual mouse pads
Drawers	16.5"W x 11.5"D x 3" and 6"H, 21" total
Drawer slides	Soft close
Laptop storage	Lockable drawer
Locks	Two zones, key or push-button
Power	Standard AC, optional battery
Drip edge	On all top and slide-out surfaces
Casters	Three locking, one tracking, 5"
Finishes	Four paint colors, two laminates





RASPBERRY MED™

PHONE

(616) 949-0080

GENERAL

info@raspberrymed.com

ORDERS

orders@raspberrymed.com

ONLINE

raspberrymed.com

Raspberry Med is a leading medical equipment manufacturer dedicated to producing high-quality medical storage solutions. Our company operates nationally, providing efficient and innovative storage products for medical facilities, laboratories, and healthcare organizations.

We take pride in offering cutting-edge storage solutions designed to meet the unique needs of our clients. At Raspberry Med, we are committed to excellence, delivering reliable products that set new standards in medical equipment manufacturing.